



Configuration FORTESSA X20 CRCT

**UV Laser
(355 nm)**




Configuration FORTESSA X20 CRCT

Violet Laser (405 nm)

↪ 300 Qdot 705 702 ↪
↪ 405 BV711 711 ↪

↪ 300 Qdot 605 603 ↪
↪ 405 BV605 603 ↪

↪ 405 BV421 421 ↪
↪ 398 Cascade Blue 420 ↪
↪ 401 Alexa Fluor 405 422 ↪
↪ 404 Pacific Blue 456 ↪
↪ 360 Dapi Viab 461 ↪



↪ 458 AmCyan 489 ↪
↪ 400 Cascade Yellow 550 ↪
↪ 400 Pacific Orange 551 ↪
↪ 405 BV510 510 ↪


↪ 300 Qdot 655 654 ↪
↪ 405 BV650 645 ↪

↪ 300 Qdot 800 792 ↪
↪ 405 BV785 785 ↪

Configuration FORTESSA X20 CRCT

Blue Laser (488 nm)

- ↻ 489 GFP 510 ↻
- ↻ 500 LysoTracker Green 510 ↻
- ↻ 493 Calcein 514 ↻
- ↻ 495 FITC 519 ↻
- ↻ 499 Alexa Fluor 488 520 ↻
- ↻ 495 CFSE 525 ↻



Configuration FORTESSA X20 CRCT

Yellow Green Laser (561 nm)

↻ 566 PE-Texas Red (ECD) (↻595) 616 ↻


↻ 351 PI 617 ↻

↻ 566 Pe-Alexa Fluor 610 (↻610) 627 ↻

↻ 566 PE-CF594 612 ↻

↻ 587 mCherry 610 ↻


↻ 566 PE-Cy5.5 (↻675) 693 ↻



Configuration FORTESSA X20 CRCT

Red Laser (640 nm)

↻650 **APC-Cy5.5** (↻675) 695↻
↻675 **Cy5.5** 695↻
↻680 **Alexa Fluor 680** 700↻
↻700 **Alexa Fluor 700** 720↻



↻631 **Alexa Fluor 633** 647↻
↻650 **APC** 660↻
↻647 **Cy5** 665↻
↻653 **Alexa Fluor 647** 668↻

↻650 **APC-Alexa Fluor 750** (↻752) 779↻
↻650 **APC-H7** 780↻
↻650 **APC-Cy7** (↻753) 785↻